

# Compella™ Bariatric bed

Advancing Bariatric Care



## Features and benefits

Maintain bariatric **patient dignity** by providing a care setting in a bed designed to resemble a standard hospital bed

Reduce caregiver strain and workload through **one-button powered width expansion** and **FlexAfoot™** allowing for quick configuration of the bed

Drive efficiency and ease of use with **SlideGuard™** designed to minimise the occurrence of the patient sliding down in bed

Ensure safe patient care with the optional **IntelliDrive™ XL Power Transport**, which offers caregivers the ability to move the bed and patient with ease

Optimise patient care with **head-of-bed alarm** and optional **3-mode bed exit alarm** that alerts caregivers of any safety concerns

The sturdy siderails, **low bed position combined with seat deflate** and ergonomic controls provide patients the confidence to safely enter and exit the bed with greater independence and dignity

The **in-bed scale** means caregivers can monitor patient weight from the comfort of the bed

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**

## Compella Bariatric bed



One-button powered width expansion



Optional: IntelliDrive XL



Point-of-Care™ caregiver controls



Optional: Patient pendant control

### Technical Specifications

Safe Working Load	500 kg
Patient weight	113-454 kg
Height adjustment	47-68 cm
Transport height	43 cm
Protection class	IPX4
<b>Bed angles</b>	
Head section	0-50°
Knee section	0-30°
Reverse Trendelenburg/Trendelenburg	+9°/-9°
<b>Sleep deck width</b>	102-127 cm
<b>Overall width</b>	
minimum (width extenders retracted)	109 cm
maximum (width extenders extended)	135 cm
<b>Sleep deck length</b>	216-234 cm
<b>Overall length</b>	
minimum (foot section retracted)	231 cm
maximum (foot section extended)	249 cm

### Surfaces

- non-powered foam
- powered air with Turn Assist
- powered air with Turn Assist and CLRT (Continuous Lateral Rotation Therapy)

<b>Castors</b>	4 x 150 mm double-band
----------------	------------------------

### Controls

standard	static patient and caregiver
optional	patient pendant

### Additional options

- 3-mode bed exit alarm
- IntelliDrive XL

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
Doc. No: 5EN126305-02, 18 May 2016

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**