



Hillrom™

Welch Allyn® ELI® Cart for Outpatient Care

Take the power of your ECG
wherever you need it.

Featuring:

1

A work surface measuring 22 in x 16 in
(55.58 cm x 40.64 cm)

2

Four swivel casters with brakes

3

An accessory bin for storage and cable
management

4

A holster for the Welch Allyn WAM™ Wireless
Acquisition Module or the Welch Allyn AM™12
Acquisition Module and patient cable



Use your ELI Cart
with the following
Hillrom ECGs:



Welch Allyn ELI 280 ECG



Welch Allyn ELI 150c ECG



Welch Allyn ELI 250c ECG

Contact your Hillrom representative to learn more.

hillrom.com

Hillrom reserves the right to make changes without notice in design, specifications and models. The only warranty Hillrom makes is the express written warranty extended on the sale or rental of its products.

Part Number: 9911-024-61

© 2019 Welch Allyn. ALL RIGHTS RESERVED. MC16438 80021081 Ver D 2019-09